

Using the Agilent xCELLigence RTCA ePacer for Functional Maturation of Human-Induced Pluripotent Stem Cell-Derived Cardiomyocytes

Authors

Xiaoyu Zhang, Jeff Li,
Miyoun Hong, and
Yama A. Abassi
Agilent Technologies, Inc.

Abstract

This application note presents a method for improving human-induced pluripotent stem cell-derived cardiomyocyte (hiPSC-CM) maturation. The Agilent xCELLigence RTCA ePacer system was used with microplates featuring planar interdigitated gold microelectrodes for directed electrical field stimulation of hiPSC-CMs. The resulting positive force-frequency relations and appropriate contractile responses to inotropic compounds show that extended electrical pacing leads to improved functional maturation of cardiomyocytes.

Introduction

One of the best-studied and characterized cell types generated from human-induced pluripotent stem cells are hiPSC-CMs, which are currently being used for drug discovery, toxicology, and investigating the underlying mechanism of cardiac disease.^{1,2} These cardiomyocytes have already been implemented as an integral *in vitro* model for toxicology and safety assessment of pharmaceutical compounds under development,^{3,4} particularly for assessment of arrhythmia.^{5,6} hiPSC-CMs display the appropriate repertoire of ion channels and currents found in human adult cardiomyocytes. Numerous studies and approaches, including patch clamp recording, field potential (FP) recording, calcium transients, impedance (IMP) measurements, and vector motion analysis have been used to validate hiPSC-CMs as a predictive model system for safety and toxicology assessment.⁷⁻¹⁰ The FDA has recently launched the Comprehensive *In Vitro* Pro-Arrhythmia Assay (CiPA) initiative to further explore and validate hiPSC-CMs for pro-arrhythmic assessment of compound risk with the aim of modifying the current regulatory guidelines for cardiac safety assessment to include hiPSC-CMs.^{11,12}

One of the main obstacles to full implementation of hiPSC-CMs for drug discovery purposes, especially those drugs that may serve to modulate the force of cardiomyocyte contraction (inotropes), has been its inherently less developed or immature phenotype.^{13,14} It is well known that the inherent structural hallmarks, calcium handling mechanism, electrophysiology, and gene expression profile of hiPSC-CMs resemble fetal rather than adult cardiomyocytes.¹⁵⁻¹⁸

From a functional perspective, hiPSC-CMs and human embryonic stem cell-derived cardiomyocytes (hESC-CMs) display a negative force-frequency relationship (NFF), whereas adult cardiomyocytes have an inherent positive force-frequency relationship (PFF).^{15,18,19} Furthermore, hiPSC-CMs display rounded and multiangular morphology, and lack the characteristic rectangular morphology of adult cardiomyocytes, which is optimized for force generation and promotion of excitation-contraction coupling.^{16,20,21} Various *in vitro* approaches, including long-term cell culture, growth factor stimulation, physical or mechanical stimulation, substrate stiffness, and electrical stimulation have been described to improve the maturation status of hiPSC-CMs.²²⁻²⁷

The maturation of cardiomyocytes can be defined from multiple aspects including structural, metabolic, electrophysiological, and functional maturation. The scope of this study is to improve hiPSC-CM maturation and evaluate it only at the functional level, specifically contractile activity of cells. To improve the maturation status of hiPSC-CMs, we have used the xCELLigence RTCA ePacer system (Figure 1A) in conjunction with microplates fabricated with planar interdigitated gold microelectrodes to perform directed electrical field stimulation of hiPSC-CMs (Figure 1B). The frequency of the applied electrical field was progressively increased, and contractility of cardiomyocytes was monitored by impedance measurement. The data demonstrate that cardiomyocytes subjected to prolonged electrical pacing display enhanced functional maturation of cardiomyocytes indicated by positive force-frequency relations and appropriate contractile responses to inotropic compounds.




Figure 1. (A) The Agilent xCELLigence RTCA ePacer is a high-throughput 6-cradle platform that consists of four components: laptop, ePacer analyzer (left), ePacer station (right), and E-Plates. (B) An image of the Agilent E-Plate CardioECR 48 and E-Plate Cardio View 96. On the E-Plate Cardio View 96, a close-up of the wells reveals the layout of the electrodes as well as the “view area”, an electrode-free gap that enables users to visually observe the cardiomyocytes in the wells.

Experimental

Materials and methods

Cell culture

hiPSC-CMs from FUJIFILM Cellular Dynamics International (FCDI/iCell cardiomyocytes²:R1017, Madison, WI, USA) were stored in liquid nitrogen until thawed and cultured according to manufacturer instructions. Each well of the Agilent E-Plate Cardio ECR 48, E-Plate Cardio 96, or E-Plate Cardio View 96 (ACEA Biosciences, Inc., a part of Agilent Technologies, San Diego, CA, USA) was coated with 50 μ L of a 1:100 diluted fibronectin (FN) solution at 10 μ g/mL (F1114, Sigma-Aldrich, St. Louis, MO, USA) and incubated at 37 °C for at least one hour. This was followed by replacing fibronectin solution with 50 μ L of prewarmed iCell cardiomyocyte-plating medium. Cells were thawed and diluted in prewarmed plating medium at 1,000,000 viable cells/mL. 50 μ L of the cell suspension was transferred using a multichannel pipette, and seeded directly onto a precoated E-Plate Cardio ECR 48 (50,000 cells/well) in a laminar hood. The plates containing iCell cardiomyocytes² (iCell CM²s) were kept in the hood at room temperature for 30 minutes, then placed and cultivated in a humidified incubator with 5% CO₂ at 37 °C. The plating medium was replaced with iCell cardiomyocyte maintenance medium for four hours postseeding. Medium change was performed every other day afterwards.

Electrical pacing of hiPSC-CMs using the ePacer

After cell seeding in the wells of the appropriate E-Plates, the E-Plates Cardio ECR 48, E-Plates Cardio 96, or E-Plates Cardio View 96 were placed in the Cardio or Cardio ECR cradles of ePacer and the impedance signals were recorded every four hours.

Following seven days of iCell CM² culture once the cells had established consistent and robust contractile activity as measured by impedance, electrical stimulation was continuously applied to the cells and progressively increased at fixed intervals. The cells were subjected to a pacing frequency starting at 1 Hz (for E-Plate Cardio 48, the stimulus setting was 0.85 V, 0.1 ms; for E-Plate Cardio 96, the stimulus setting was 1 V, 0.11 ms), followed by 1.5 Hz (for E-Plate Cardio ECR 48, the stimulus setting was 0.9 V, 0.12 ms; for E-Plate Cardio 96, the stimulus setting was 1.1 V, 0.12 ms) and ended at 2 Hz (for E-Plate Cardio ECR 48, 0.95 V, 0.18 ms; for E-Plate Cardio 96, the stimulus setting was 1.1 V, 0.15 V). The duration of each applied pacing frequency was five days. The stimulus setting mainly depends on cell status and properties, which could vary. The entire pacing process takes approximately two weeks. Alternatively, a one-week maturation protocol can be used seven days postseeding. Cells are subjected to electrical stimulation at 1 Hz for one day, followed by a consecutive three-day pacing at 2 Hz using a similar stimulus setting to the two-week pacing protocol.

Chemical reagents

All the chemical reagents were purchased from Tocris (Minneapolis, MN, USA) and Sigma-Aldrich (St. Louis, MO, USA). 1,000-fold chemical stock solutions were prepared in DMSO and stored at -20 °C. The serial diluted chemicals (1,000-fold) were further prepared in DMSO immediately before compound addition. The 10-fold final dilution of the chemicals was prepared with culture medium for single use only. The final concentration of DMSO in the treated well was 0.1%.

The RTCA ePacer and assay principle

The xCELLigence RTCA ePacer is a dual-mode instrument that includes both real-time monitoring of hiPSC-CM viability and contraction as well as directed electrical pacing of hiPSC-CM. It consists of a workstation (computer and software), an electronic analyzer, a plate station that is placed inside a CO₂ incubator, and E-Plates, which are specialized electronic microplates (Figure 1B) placed within the cradles of the ePacer. Embedded within the bottom of the E-Plate wells are interdigitated gold microelectrodes, which noninvasively monitor cell impedance signal, providing both a measure of hiPSC-CM viability and, under high-frequency data acquisition mode, a measure of contractility (Figure 2B). The cellular impedance signal is recorded at a user-defined time interval (minutes and hours), and is reported using a unitless parameter called Cell Index. During the electrical pacing, the electrical pulses are directly applied to the cells through all the IMP electrodes on the bottom of the wells. For most cardiomyocytes, the length of each electrical pulse used by IMP electrode is less than 4 ms, which allows the contractile activities of cells to be immediately captured and recorded while the cells are being paced by IMP electrodes. The exact conditions for electrical pacing are dependent on the cell type, the inherent beating frequency, and the experimental context.

Three types of E-plates are compatible with the cradles of ePacer, including E-Plate Cardio ECR 48 (48-well plate), E-Plate Cardio 96, and E-Plate Cardio View 96 (96-well plate). The E-Plate Cardio View 96 contains a "view area" allowing imaging of the cells (Figure 1B).

The ePacer was developed and designed to be an integral component of a workflow that can be used in conjunction with both the xCELLigence Cardio and CardioECR systems as well as other orthogonal readouts such as calcium imaging, high-content imaging, fluorescent voltage dye imaging, and plate reader assays. Upon completion of the electrical pacing regimen on the ePacer, the plates containing the cells can either be transferred to xCELLigence Cardio or CardioECR instruments and continue to be monitored and treated with compounds. Alternatively, the cells within the plates can be processed appropriately and transferred to other types of instruments, as outlined in Figure 2C.





Figure 2. (A) An image of the Agilent xCELLigence E-Plate Cardio 96, which shows the interdigitated impedance (IMP) electrode array at the bottom of the wells. (B) Cardiomyocyte contraction and relaxation induce changes in cell morphology and cell attachment, which can be detected by impedance readout. (C) Full assay workflow using the Agilent xCELLigence RTCA ePacer. (1) Simply pace the cells on the ePacer, (2) assess maturity using the ePacer software, and (3) integrate into an assay platform of your choice.

Results and discussion

Using the ePacer for assessment of pacing efficiency and cell quality during chronic electrical stimulation

Figure 3A shows the dynamic changes in attachment and viability of iCell CM²s after seeding in E-Plate Cardio 96 for seven days, as measured by impedance and displayed as an overall Cell Index curve. Each data point in the overall Cell Index curve contains a recording of 15 seconds of high frequency impedance data, which represents the actual contraction of hiPSC-CMs generated at different time points along the cell culture. Real-time monitoring of both overall Cell Index and contraction in terms of rate and amplitude is a very useful parameter for quality assessment of hiPSC-CM. As the cells attach, they form a monolayer that is both physically and electrically coupled, culminating in synchronous and stable beating of the cells within seven days after cell seeding. To demonstrate the precision of directed electrical pacing, iCell CM²s were seeded in different wells and were subjected to electrical pacing at increasing frequencies. As shown in Figure 3B, the iCell CM²s were able to precisely follow the increasing frequency of the directed electrical pacing by displaying increased beating rate. Figure 3C shows a direct linear relationship between the pacing frequency and observation of increased beating rate.

With respect to directed electrical pacing, there are two important parameters to achieve optimal functional cells. These parameters include pacing efficiency as well as overall viability of hiPSC-CMs in response to directed electrical pacing. It is critical that the hiPSC-CMs follow the input pacing stimulus in a precise and uninterrupted manner. The ePacer automatically calculates and displays the number of wells that follow the input pacing




Figure 3. (A) The overall Cell Index curves were recorded in real time after cell seeding. Each data point in the Cell Index curve contains a recording of 15 seconds of IMP waveform data, which reflects the actual contraction of the cardiomyocytes, as shown in the figure. (B) These impedance waveforms show the cell beating rate response to the pacing workflow of the Agilent xCELLigence RTCA ePacer, from spontaneous beating to the 2 Hz target pacing frequency. (C) A linear relationship is observed between the pacing frequency and the cell beating rate. Continued next page.

stimulus (pacing efficiency, displayed as a percentage). This real-time display of pacing efficiency is important for making appropriate adjustments to input stimulus (Figure 3D). The data collected from two different test sites show high pacing efficiency in response to different input pacing stimuli (Table 1). This demonstrates that the ePacer is

optimally designed to achieve consistent and reliable beating rates, which are important for functional maturation of hiPSC-CMs.

Cell quality and viability are other key factors affecting cell responses to the treatment. Typically, chronic electrical pacing should be performed under conditions that are not stressful to

the cell and do not adversely impact cell viability. The ePacer monitors and displays overall cell health in real time both under pacing or nonpacing conditions. Figure 3E shows that there was no significant decrease in overall Cell Index obtained from paced cells compared to the nonpaced cells.





Figure 3. (D) The top section of this screenshot from the ePacer software shows the pacing efficiency of the cells when the electrical pacing is applied. The bottom section of the screenshot shows the cells contractile activity of each well in the plate. (E) The overall Cell Index was measured in nonpaced (CTRL) cells and electrically paced cells.

Table 1. Summary of the pacing capture efficiency of the Agilent xCELLigence RTCA ePacer during 2-week electrical pacing from two different test sites. The beating rate at different stages of the pacing process was presented as mean \pm STDEV (n = 48 per site).

| Beating Rate | Day 5_1Hz (60 beats) | | Day 3/5_1.5 Hz (90 beats) | | Day 5_2 Hz (120 beats) | |
|--------------|----------------------|-------------------|---------------------------|-------------------|------------------------|-------------------|
| | Mean \pm Std. Dev. | Pacing Efficiency | Mean \pm Std. Dev. | Pacing Efficiency | Mean \pm Std. Dev. | Pacing Efficiency |
| Site 1 | 59.8 \pm 0.7 | 100% | 89.9 \pm 0.1 | 100% | 118.1 \pm 3.7 | 90% |
| Site 2 | 59.9 \pm 0.12 | 100% | 90.6 \pm 1.7 | 100% | 121.4 \pm 1.9 | 100% |

Directed electrical pacing reverses the BAmplitude and beating rate relationship and improves the functional response to inotropic compounds

One of the hallmarks of hiPSC-CM immaturity is the inherent negative force-frequency relationship.¹⁸ To determine whether long-term electrical stimulation through the planar IMP electrodes impacts the functional maturation of hiPSC-CMs, we sought to determine the FF relationship in both nonpaced and electrically paced cells. Since measured impedance amplitude (BAmplitude) can be used as a surrogate for force of contraction, we evaluated the BAmplitude and beating rate relationship in the cells as a surrogate for the FF relationship. We cultured hiPSC-CMs in the E-Plates, and subjected half of the plate to directed electrical pacing (paced group) while maintaining the remaining half of the E-Plate without electrical stimulation (control group). The paced group was subjected to chronic electrical pacing, increasing the pacing frequency from 1 Hz (week 1) to 1.5 Hz (week 2) and finally 2 Hz (week 3) (Figure 4A). At the end of each pacing regimen, the hiPSC-CMs were subjected to acute electrical pacing with increasing

frequencies (0.75, 1, 1.5, 2) as described previously. The data show that before applying stimulation to the cells after 1 week in culture (week 0), both paced and nonpaced cells displayed a negative BAmplitude-frequency relationship (Figure 4B and C). Continued stimulation at weeks 1 and 2 seemed to decrease the slope of the BAmplitude-frequency relationship, and at the end of week 3, a positive BAmplitude-frequency relationship was observed in the electrically paced group (Figure 4C). As a control, nonpaced cells were subjected to acute electrical stimulation, and continued to show negative BAmplitude-frequency relationship (Figure 4B). A shorter pacing protocol was also tested as described in the materials and methods section; it also showed a positive BAmplitude-frequency relationship (Figure 4D)

To determine if chronic electrical pacing of hiPSC-CM, as previously described, improved functional response of the cells to inotropic compounds, we first tested isoproterenol (ISO), a well known positive inotropic and chronotropic compound. As shown in Figure 5A, ISO (100 nM) treatment resulted in increased beating rate for both paced and nonpaced cells.

However, the inotropic effect of ISO could only be observed in electrically paced cells, as demonstrated by an increase in the BAmplitude (52% \pm 14% SD) (shown in Figure 5B), while nonpaced cells showed a slight decrease in BAmplitude, as shown previously. We further expanded the functional test in paced cells by testing more positive inotropic compounds, with different mechanisms of action, as well as a negative inotropic compound. Table 2 shows that positive inotropic compounds profoundly increased beating amplitude only in paced cells and not in nonpaced cells (data not shown). Conversely, a negative inotropic compound, isradipine, decreased beating amplitude in paced and nonpaced cells. In summary, chronic electrical pacing converts hiPSC-CM response from an inherent negative BAmplitude-frequency relationship to a positive BAmplitude-frequency relationship and improves the cell contractile responses to positive inotropes, suggesting that long-term electrical pacing improves the maturation status of hiPSC-CMs at the functional level.




Figure 4. The timeline of chronic electrical pacing protocol (A). Beating amplitude (BAmp) and beating rate relationship was obtained from (B) nonpaced cells and (C) paced cells in week 0, the time before chronological pacing started; week 1, the time before 1 Hz pacing ended; week 2, the time before 1.5 Hz pacing ended; and week 3, the time before 2 Hz pacing ended.

Conclusion

hiPSC-CMs are derived from reprogrammed somatic cells. Numerous lines of evidence, including morphological, structural, gene and protein expression, and functional, indicate that these cells display the hallmarks of fetal rather than adult cardiomyocytes²⁹, and various *in vitro* approaches have failed to recapitulate the natural cardiomyocytes development program. One of the key functional consequences of the inherent immaturity and lack of sophisticated myocyte structure of hiPSC-CM is generation of a negative force-frequency (NFF) relationship.^{18,19} Using IMP measurement as a surrogate for contractility and directed electrical pacing to control the rate of contraction, our results demonstrate that hiPSC-CMs display a negative BAmplitude-frequency relationship (Figure 4A), confirming previous observations that hiPSC-CM display a NFF relationship.

To improve the contractile maturation of hiPSC-CMs, we developed a novel instrument, the Agilent xCELLigence RTCA ePacer, which uses microplates embedded with electrodes to perform directed electrical pacing of the cells in a highly controlled and consistent fashion for both acute and chronic stimulation. In this study, we used the ePacer to perform chronic electrical pacing of hiPSC-CMs, ranging from 1 to 2 Hz. Our results confirm the earlier findings that electrical pacing of hiPSC-CM can lead to further maturation.^{27,30} In addition, our data clearly demonstrate for the first time that the conditions we used for electrical pacing reverse the

negative BAmplitude-frequency relationship of hiPSC-CM, and result in a positive BAmplitude-frequency relationship, indicating that electrical stimulation can lead to further functional maturation of hiPSC-CM (Figure 4). Furthermore, the

mature FF relationship observed in chronically paced cells allows the cells to respond to positive inotropic compounds in an appropriate and dose-dependent manner (Figure 5B and Table 2).




Figure 5. Isoproterenol, a known positive inotrope, was added to both nonpaced and paced cells. (A) The % change in beating rate 30 minutes after isoproterenol (ISO) addition to nonpaced beating cells (blue bar) and chronically paced cells (orange bar) at different concentrations. (B) The % change in BAmplitude obtained from nonpaced beating cells (blue bar) and chronically paced cells (orange bar). The data were represented by mean \pm STDEV ($n \geq 3$).

In addition to the powerful electrical pacing function, the ePacer system comes equipped with its own software, allowing both real-time monitoring of the cells in terms of viability and contractility, and allowing precise programming of all the pacing parameters in terms of magnitude, duration, and shape of the input stimulus. It also allows defined and controlled pacing protocols while providing the throughput and flexibility needed to screen compounds or multiple replicates. The automatic FF relationship evaluation feature of the ePacer software helps the user to determine when the paced cells have developed a positive force frequency relationship, which is an attribute of matured cardiomyocytes.

The compatibility of the E-Plate Cardio View 96 with orthogonal readouts such as calcium and high-content imaging allows integration of this system into various workflows depending on the experimental needs.

Table 2. Summary of contractile responses of chronic electrically paced iCell CM²s. After long-term pacing, the cells were exposed to both positive and negative inotropic compounds. The % change of BAmp after compound addition was further calibrated to the time-matched DMSO control wells. The data were represented by mean \pm STDEV (n \geq 3).

| | Bay K8644 (67 nM) | Pimobendan (10 μM) | Milrinone (30 μM) | Omecamtiv mecarbil (600 nM) | Digoxin (100 nM) | Isradipine (10 nM) |
|-------------------------|---|--|---|--|------------------------------|---|
| MOA | L-type Ca ²⁺ ion channel activator | Ca ²⁺ sensitizer | Phosphodiesterase-3 inhibitor | Myosin activator | Na-K ATP exchanger inhibitor | L-type Ca ²⁺ ion channel inhibitor |
| Inotropic Effect | Positive | Positive | Positive | Positive | Positive | Negative |
| % Change (BAmp) | 32 \pm 5% | 31 \pm 6% | 30 \pm 12% | 58 \pm 15% | 175 \pm 35% | -47 \pm 7% |

References

1. Duelen, R.; Sampaolesi, M. Stem Cell Technology in Cardiac Regeneration: A Pluripotent Stem Cell Promise. *EBioMedicine* **2017**, *16*, 30–40.
2. Iglesias-Garcia, O.; Pelacho, B.; Prosper, F. Induced Pluripotent Stem Cells as a New Strategy for Cardiac Regeneration and Disease Modeling. *J. Mol. Cell Cardiol.* **2013**, *62*, 43–50.
3. Takasuna, K. *et al.* Comprehensive *in vitro* Cardiac Safety Assessment Using Human Stem Cell Technology: Overview of CSAHi HEART Initiative. *J. Pharmacol. Toxicol. Methods* **2017**, *83*, 42–54.
4. Gintant, G.; Sager, P. T.; Stockbridge, N. Evolution of Strategies to Improve Preclinical Cardiac Safety Testing. *Nat. Rev. Drug Discov.* **2016**, *15*, 457–71.
5. Ando, H. *et al.* A New Paradigm for Drug-Induced Torsadogenic Risk Assessment Using Human iPSC Cell-Derived Cardiomyocytes. *J. Pharmacol. Toxicol. Methods* **2016**, *84*, 111–127.
6. Guo, L. *et al.* Estimating the Risk of Drug-Induced Proarrhythmia Using Human Induced Pluripotent Stem Cell-Derived Cardiomyocytes. *Toxicol. Sci.* **2011**, *123*, 281–9.
7. Qu, Y.; Vargas, H. M. Proarrhythmia Risk Assessment in Human Induced Pluripotent Stem Cell-Derived Cardiomyocytes Using the Maestro MEA Platform. *Toxicol. Sci.* **2015**, *147*, 286–95.
8. Zhang, X. *et al.* MultiParametric Assessment of Cardiomyocyte Excitation-Contraction Coupling Using Impedance and Field Potential Recording: A Tool for Cardiac Safety Assessment. *J. Pharmacol. Toxicol. Methods* **2016**, *81*, 201–16.
9. Guo, L. *et al.* Refining the Human iPSC-Cardiomyocyte Arrhythmic Risk Assessment Model. *Toxicol. Sci.* **2013**, *136*, 581–94.
10. Hayakawa, T. *et al.* Image-Based Evaluation of Contraction-Relaxation Kinetics of Human-Induced Pluripotent Stem Cell-Derived Cardiomyocytes: Correlation and Complementarity with Extracellular Electrophysiology. *J. Mol. Cell Cardiol.* **2014**, *77*, 178–91.
11. Gintant, G. *et al.* The Evolving Roles of Human iPSC-Derived Cardiomyocytes in Drug Safety and Discovery. *Cell Stem Cell* **2017**, *21*, 14–17.
12. Fermini, B. *et al.* A New Perspective in the Field of Cardiac Safety Testing through the Comprehensive In Vitro Proarrhythmia Assay Paradigm. *J. Biomol. Screen.* **2015**, *21*, 1–11.
13. Ravenscroft, S. M. *et al.* Cardiac Non-Myocyte Cells Show Enhanced Pharmacological Function Suggestive of Contractile Maturity in Stem Cell Derived Cardiomyocyte Microtissues. *Toxicol. Sci.* **2016**, *152*, 99–112.
14. Pillekamp, F. *et al.* Contractile Properties of Early Human Embryonic Stem Cell-Derived Cardiomyocytes: Beta-Adrenergic Stimulation Induces Positive Chronotropy and Lusitropy But Not Inotropy. *Stem Cells Dev.* **2012**, *21*, 2111–21.
15. Binah, O. *et al.* Functional and Developmental Properties of Human Embryonic Stem Cells-Derived Cardiomyocytes. *J. Electrocardiol.* **2007**, *40*, S192–6.
16. Hartman, M. E.; Dai, D. F.; Laflamme, M. A. Human Pluripotent Stem Cells: Prospects and Challenges as a Source of Cardiomyocytes for *in vitro* Modeling and Cell-Based Cardiac Repair. *Adv. Drug Deliv. Rev.* **2015**, *96*, 3–17.
17. Robertson, C.; Tran, D. D.; George, S. C. Concise Review: Maturation Phases of Human Pluripotent Stem Cell-Derived Cardiomyocytes. *Stem Cells* **2013**, *31*, 829–37.
18. Germanguz, I. *et al.* Molecular Characterization and Functional Properties of Cardiomyocytes Derived from Human Inducible Pluripotent Stem Cells. *J. Cell Mol. Med.* **2011**, *15*, 38–51.
19. Dolnikov, K. *et al.* Functional Properties of Human Embryonic Stem Cell-Derived Cardiomyocytes: Intracellular Ca²⁺ Handling and the Role of Sarcoplasmic Reticulum in the Contraction. *Stem Cells* **2006**, *24*, 236–45.
20. Zhang, J. *et al.* Functional Cardiomyocytes Derived from Human Induced Pluripotent Stem Cells. *Circ. Res.* **2009**, *104*, e30–41.
21. Itzhaki, I. *et al.* Calcium Handling in Human Induced Pluripotent Stem Cell Derived Cardiomyocytes. *PLoS One* **2011**, *6*, e18037.
22. Lundy, S. D. *et al.* Structural and Functional Maturation of Cardiomyocytes Derived from Human Pluripotent Stem Cells. *Stem Cells Dev.* **2013**, *22*, 1991–2002.
23. Yang, X. *et al.* Tri-Iodo-L-Thyronine Promotes the Maturation of Human Cardiomyocytes-Derived from Induced Pluripotent Stem Cells. *J. Mol. Cell Cardiol.* **2014**, *72*, 296–304.
24. Herron, T. J. *et al.* Extracellular Matrix-Mediated Maturation of Human Pluripotent Stem Cell-Derived Cardiac Monolayer Structure and Electrophysiological Function. *Circ. Arrhythm. Electrophysiol.* **2016**, *9*, e003638.

25. Aigha, I.; Raynaud, C. Maturation of Pluripotent Stem Cell Derived Cardiomyocytes: The New Challenge. *Glob. Cardiol. Sci. Pract.* **2016**, e201606.
26. Ruan, J. L. *et al.* Mechanical Stress Promotes Maturation of Human Myocardium from Pluripotent Stem Cell-Derived Progenitors. *Stem Cells* **2015**, *33*, 2148–57.
27. Chan, Y. C. *et al.* Electrical Stimulation Promotes Maturation of Cardiomyocytes Derived from Human Embryonic Stem Cells. *J. Cardiovasc. Transl. Res.* **2013**, *6*, 989–99.
28. Joulin, O. *et al.* Cardiac Force-Frequency Relationship and Frequency-Dependent Acceleration of Relaxation are Impaired in LPS-Treated Rats. *Crit. Care* **2009**, *13*, R14.
29. Yang, X.; Pabon, L.; Murry, C. E. Engineering Adolescence: Maturation of Human Pluripotent Stem Cell-Derived Cardiomyocytes. *Circ. Res.* **2014**, *114*, 511–23.
30. Nunes, S. S. *et al.* Biowire: A Platform for Maturation of Human Pluripotent Stem Cell-Derived Cardiomyocytes. *Nat. Methods* **2013**, *10*, 781-7.
31. Scott, C. W. *et al.* An Impedance-Based Cellular Assay Using Human iPSC-Derived Cardiomyocytes to Quantify Modulators of Cardiac Contractility. *Toxicol. Sci.* **2014**, *142*(2), 331–8.

www.agilent.com/chem/ePacer

For Research Use Only. Not for use in diagnostic procedures.

This information is subject to change without notice.

© Agilent Technologies, Inc. 2019
Printed in the USA, December 6, 2019
5994-1552EN

Contact OLS OMNI Life Science - Your Partner in Cell Research

www.ols-bio.com

