

CASY^{VIVO} Cell Counter & Analyzer - Powering Breakthroughs in Cell Research

Standardizing Biocompatibility Assays via Precise Cell Counting

Woitschach et al. (2021). Innate Immune Cell Interaction with polymer-coated Silicone; Materials, 14(10):2474. DOI: 10.3390/ma14102474.

Immunology; Monocytes; Biocompatibility	
Index	IM12
Standardization	X
Counting	X
Viability	
Volume	

The Challenge:

Assessing the biocompatibility of new polymer-coated silicone (LSR-PMPC, LSR-PUR) by quantifying innate immune cell (monocyte/macrophage) adhesion and activation on these surfaces.

CASY's Contribution:

CASY determined the precise cell number and viability of isolated human monocytes. This standardized the cell seeding density (1×10^5 cells/well) for all adhesion and cytokine assays. It was also used to count the adherent T-cells, proving that LSR-PMPC surfaces significantly reduced monocyte adhesion⁴⁴⁴

Key Benefits to Researchers:

- Standardization:** Provided accurate, label-free counting of primary monocytes post-isolation, ensuring a precise and consistent T-cell seeding density for all experiments.
- Accuracy:** Enabled quantitative, objective measurement of cell adhesion by accurately counting the number of adherent T-cells remaining on each material after 2 hours.
- Validation:** This precise quantification was essential for validating the key finding: that the zwitterionic (LSR-PMPC) and polyurethane (LSR-PUR) coatings significantly reduce monocyte adhesion, indicating improved biocompatibility.

CASY Analysis of Monocyte Adhesion on TPU and LSR Surfaces




Figure 2. Decreased monocyte adhesion on thermoplastic polyurethane (TPU) and liquid silicone rubber (LSR) with only minor effects on cell viability and metabolic activity. Then, 5×10^5 cells were incubated on the materials ... **(C)** Non-attached cells were determined after 12 h of incubation by collecting the supernatants and counting the non-attached cells with Casy cell counter.