

# A New Way to Monitor Virus-Mediated Cytopathogenicity

## Author

Brandon Lamarche,  
Joyce Velez, and Leyna Zhao  
Agilent Technologies, Inc.

## Introduction

One of the most important procedures in virology is the measurement of viral cytopathic effects (CPEs). The plaque assay has long been the gold standard for quantifying CPEs by providing a direct readout of the number or concentration of infectious viral particles in a sample. In this technique, a confluent monolayer of host cells is infected with varying dilutions of the virus and is overlaid with a semisolid material, such as dilute agarose gel. When an infected cell lyses, the overlay material prevents the released virions from diffusing through the medium and infecting distal sites. However, progeny virions can gain access to neighboring cells in the immediate vicinity. In this manner, infection and lysis spread laterally in two dimensions and produce a cell-free plaque in the middle of an otherwise confluent group of cells. Depending on the host cell type and the virus, accurate recognition and counting of plaques may require staining cells with a dye, such as crystal violet.

A very low multiplicity of infection (MOI) can be achieved using serial dilutions of virus. Under this condition, each cell that gets infected will be infected by just one virion. After counting the number of plaques in a well, and accounting for the dilution factor used in preparing the virus inoculum, the concentration of virus (titer) in the original sample can be calculated. Titers are usually reported as the number of plaque forming units (PFU) per unit of volume.

Depending on the virus and host cells being studied, viral plaque formation can take anywhere from days to weeks to be detectable. A single endpoint plaque assay provides no information about the onset of CPE or the kinetics of virus-mediated cytotoxicity. Different cell types and cell densities, as well as viral strains, serotypes, and mutations can cause plaque formation rates and sizes to vary dramatically. Thus, the suboptimal selection of a single assay endpoint can result in inaccurate calculation of viral titer and lytic activity. The definition and manual counting of plaques by visual inspection can also be highly subjective.

Recent peer-reviewed studies of oncolytic viruses (Dyer *et al.*, 2017<sup>2</sup> and Fajardo *et al.*, 2017<sup>5</sup>) and cancer vaccines (Cross *et al.*, 2015<sup>1</sup> and Pham *et al.*, 2014<sup>6</sup>) have demonstrated that the Agilent xCELLigence real-time cell analysis (RTCA) system is a powerful tool for evaluating both virus concentration and cytotoxicity kinetics. It uses a simple, fast, and reproducible workflow. Microelectronic biosensors embedded in the bottom of microplate wells enable the RTCA assay to offer dynamic, real-time, label-free, and noninvasive analysis of cellular events, such as virus-mediated cytolysis. The progeny viruses released from a lysed cell are free to diffuse through the media and infect distant target cells because an agarose overlay is not used in the RTCA procedure. This unhindered spread of virus throughout the entire well results in the rapid lysis of all cells, providing a quantification of viral titer much more quickly than a plaque assay.

xCELLigence impedance measurements are automatically recorded at a user-defined frequency and are plotted by the xCELLigence software using the dimensionless parameter known as Cell Index (CI).

Key benefits of the xCELLigence RTCA systems for monitoring virus-mediated cytopathogenicity:

- **Label free:** No dyes required.
- **Fast:** Read an entire 96-well plate in less than 10 seconds.
- **Real time:** Quantitative monitoring of both fast (hours) and slow (days to weeks) CPE.
- **Easy workflow:** No gel pouring. Requires only the addition of virus to host cells.
- Accurate, precise, and highly reproducible.

- **Automatic data plotting:** The intuitive xCELLigence software enables easy data display and objective analysis, precluding the subjective data vetting that is common to plaque assays.

This application note describes the experimental setup for assessing vesicular stomatitis virus (VSV)-mediated cytotoxicity of Vero E6 cells and HEK 293 cells using an xCELLigence instrument. The protocol shows the identification of cell proliferation kinetics as well as the optimal time point for viral infection with different cell seeding densities. The assay overcomes many of the limitations of single-point plaque assays, and provides direct evidence that RTCA can offer a comprehensive and reliable evaluation of viral cytopathogenicity.

## Materials and methods

### Cells

Cells were cultured in a standard humidified incubator at 37 °C with 5% CO<sub>2</sub> saturation. Vero E6 (obtained from the ATCC) is an African green monkey kidney-derived cell line with deficiency of the type I interferon genes. HEK 293 (obtained from Microbix Biosystems) is a human embryonic kidney cell line with an intact interferon system. Both adherent cell lines were maintained in Dulbecco's Modified Eagle Medium (DMEM), supplemented with 10% fetal calf serum, 2 mM glutamine, and 1% penicillin/streptomycin.

### Virus

The VSV, Serotype Indiana, was grown and titrated on Vero E6 cells at 37 °C with 5% CO<sub>2</sub>.

### Cell proliferation assays

For real-time cell analysis, 100 µL of growth media was added to each well of the Agilent E-Plate 96 to obtain background readings. For each cell type, a sequential 1:1 dilution series, with seven different cell numbers ranging from 50,000 to 781 cells/well, were resuspended in 100 µL of media, then seeded into the E-Plate 96. The E-Plates containing cells were incubated for 30 minutes at room temperature, and placed on the RTCA single plate (SP) station, located in the cell culture incubator. Cell attachment, spreading, and proliferation were monitored every 30 minutes using the RTCA SP. Measured impedance recordings from cells in each individual well on the E-Plate 96 were automatically converted to CI values by the RTCA software.


### Assessment of virus-mediated cytopathogenicity

For viral studies, 25,000 cells/well and 12,500 cells/well of each cell line were seeded into each well of an E-Plate 96. When the cells reached confluency (25,000 cells/well) or were still in the growth phase (12,500 cells/well), after 20.5 hours for Vero E6 cells and 68.5 hours for HEK 293 cells, they were infected by VSV. This was done by removing the E-Plate 96 from the RTCA SP station and adding 800,000 ("high MOI") or 80,000 ("low MOI") resuspended in 10 µL of growth media to the wells. As the control, eight wells were mock-infected by adding 10 µL growth media only. The E-Plate 96 was then placed immediately back into the RTCA SP station in the incubator and the CI values were measured every 15 minutes for up to 190 hours.

## Results and discussion

### Dymanic monitoring of cell proliferation

To identify the optimal time point for viral infection, a cell proliferation analysis was performed with Vero E6 and HEK 293 cells. Suitable time points for virus infection were defined at 20.5 hours for Vero E6 cells and 68.5 hours for HEK 293 cells (Figures 1A and 1B). At these time points, cells were either in the growth phase when 12,500 cells had been used for seeding, or in the early stationary phase when 25,000 cells had been used for seeding. Therefore, viral cytopathogenicity was monitored in either the growth phase or the early stationary phase.




**Figure 1.** Dynamic monitoring of cell proliferation. Cells were seeded in the E-Plate 96 and continuously monitored by measuring CI to identify a suitable time point for addition of virus (growth or early stationary phase). The adhesion, spreading, and proliferation of (A) Vero E6 cells and (B) HEK 293 cells were dynamically monitored every 30 minutes using the RTCA SP instrument. Colored curves indicate the different cell numbers seeded per well in an E-Plate 96 (from left to right): red, 50,000; green, 25,000; blue, 12,500; magenta, 6,250; cyan, 3,125; orange, 1,562; dark green, 781; olive green, medium control (without cells).

### VSV cytopathogenicity profile using Vero E6 cells

At 20.5 hours after seeding (based on the dynamic monitoring of cell proliferation), Vero E6 cells either in growth phase or early stationary phase were infected with VSV using two different MOIs.

When Vero E6 cells were infected with VSV during the growth phase, there was a clear correlation between the amount of virus used for infection and the onset of the virus-mediated CPE (Figure 2). After infection with a low MOI (80,000 PFU VSV), the cells continued to grow for 15 hours (Figure 2A, blue curve), similar to mock-infected cells (Figure 2A, green curve). The CI values then decreased, indicating that the cells were dying as a consequence of VSV replication. In contrast, mock-infected cells continued to grow. At 24 hours after infection, the CI values had decreased to 50% of the maximum value (CI50), then continued to decline to zero, indicating complete cell death in the infected culture. In contrast, the CI of Vero E6 cells infected with a high MOI (800,000 PFU VSV) started to decline at 4 hours post infection (Figure 2A, red curve), and the CI50 was reached after 11 hours.

Very similar results were obtained when confluent Vero E6 cells were infected (Figure 2B). The CI50 was reached at 10 hours postinfection (high MOI, Figure 2B, red curve) and 19 hours postinfection (low MOI, Figure 2B, blue curve), respectively. Again, complete death of the infected cultures was observed, as indicated by the decrease of CI values to zero.



**Figure 2.** Dynamic monitoring of Vero E6 cells during VSV infection. (A) Normalized Cell Index values of growing cells, (B) Normalized Cell Index values of confluent cells. The virus-mediated effect on adhesion, spreading, and proliferation of the cells was monitored by measuring cell impedance every 15 minutes using the RTCA SP instrument. Time of addition of virus at 20.5 hours is indicated by the black vertical line. The time point when the CI value had decreased to 50% of the maximum (CI50) value is indicated by the dotted orange lines. Green curves: control (no viral infection); blue curves: 80,000 PFU VSV; red curves: 800,000 PFU VSV.

### VSV cytopathogenicity profile using HEK 293 cells

Based on the dynamic monitoring of cell proliferation, HEK 293 cells in the growth phase and in early stationary phase were infected with VSV 68.5 hours after seeding (Figures 3A and 3B).

HEK 293 cells showed a different response compared to Vero E6 cells when infected with VSV. HEK 293 cells in the growing phase were much more sensitive to VSV infection. This is indicated by the drop in CI values to the CI50 value 6 hours after infection when a high MOI was used (Figure 3A, red curve). Cells infected with a low MOI reached the CI50 12 hours after infection (Figure 3A, blue curve). A completely different result was obtained when confluent HEK 293 cells were infected (see Figure 3B). Confluent cells infected with a high MOI exhibited a drop in CI values, similar to growing cells (Figure 3B, red curve). Cells infected with a low MOI appear to be completely resistant to VSV infection, exhibiting CI values virtually identical to mock-infected cells (Figure 3B, blue and green curves). Considering the different responses to the VSV infection, the main difference between Vero E6 and HEK 293 cells was the ability to produce type I interferons. Vero E6 cells are devoid of the interferon genes.<sup>3</sup>


As a consequence, Vero E6 cells could not upregulate the expression of interferon-induced antiviral active proteins, such as MxA and OAS/RNaseL, in response to viral infections.

In contrast, HEK 293 cells possess an intact interferon system. During viral infection, they produced interferons that activated the JAK/STAT signaling pathway in an autocrine and paracrine manner. As a result, the expression of antiviral active proteins was initiated and an antiviral state was established. The observed resistance of confluent

HEK 293 cells to VSV infection with a low MOI may have been due to the antiviral response mounted by their interferon system.

Confluent cells may also represent a suboptimal environment for VSV replication because they have reduced metabolic activity compared to growing cells. In line with this hypothesis was the observation that growing HEK 293 cells are much more sensitive to VSV, independent of the MOI used for infection.

However, the VSV-M protein has been known to counteract the interferon system by inhibiting host RNA and protein synthesis. This contributed to the shutoff of host-directed gene expression.<sup>4</sup> Therefore, the observed differences in response to VSV infection in growing or confluent Vero E6 and HEK 293 cells and the dependency of the outcome on cell number (and VSV MOI) most likely reflected the interplay of cellular antiviral mechanisms and viral countermeasures.



**Figure 3.** Dynamic monitoring of HEK 293 cells during viral infection. (A) Normalized Cell Index values of growing cells, (B) Normalized Cell Index values of confluent cells. The virus-mediated effect on adhesion, spreading, and proliferation of the cells was dynamically monitored every 15 minutes using the RTCA SP instrument. Time of addition of virus at 68.5 hours is indicated by the black vertical line. The time point when the CI had decreased to 50% of the maximum value (CI50) is indicated by the dotted orange lines. Green curves show the control (no viral infection); blue curves: 80,000 PFU VSV; red curves: 800,000 PFU VSV.

In contrast to conventional endpoint assays, real-time cell analysis using the xCELLigence RTCA system offers continuous monitoring of virus-host interactions to better define the responses where viral or cellular activities are more dominant.

## References

1. Cross, *et al.* Therapeutic DNA Vaccination Against Colorectal Cancer by Targeting the MYB Oncoprotein. *Clin. Transl. Immunology* **2015**, 4(1), e30.
2. Dyer, *et al.* Oncolytic Group B Adenovirus Enadenotucirev Mediates Non-apoptotic Cell Death with Membrane Disruption and Release of Inflammatory Mediators. *Molecular Therapy: Oncolytics* **2017**, 4, 18–30.
3. Emeny, J. M.; Morgan, M. J. Regulation of the Interferon System: Evidence That Vero Cells Have a Genetic Defect in Interferon production. *J. Gen. Virol.* **1979**, 43(1), 247–252.
4. Ferran, M. C.; Lucas-Lenard, J. M. The Vesicular Stomatitis Virus Matrix Protein Inhibits Transcription From the Human Beta Interferon Promoter. *J. Virol.* **1997**, 71(1), 371–377
5. Fajardo, *et al.* Oncolytic Adenoviral Delivery of an EGFR-Targeting T Cell Engager Improves Antitumor Efficacy, **2017**. doi: 10.1158/0008-5472.CAN-16-1708
6. Pham, *et al.* A simple *in vitro* Method for Evaluating Dendritic Cell-Based Vaccinations. *Onco Targets Ther.* **2014**, 7, 1455–64.

For Research Use Only. Not for use in diagnostic procedures.

RA.5126273148  
DE.3109143519

This information is subject to change without notice.

© Agilent Technologies, Inc. 2020, 2021  
Printed in the USA, February 2, 2021  
5994-1764EN

## Contact OLS OMNI Life Science - Your Partner in Cell Research

OLS OMNI Life Science GmbH & Co. KG  
Karl-Ferdinand-Braun-Straße 2  
28359 Bremen, Germany  
Phone: +49 421 27 61 69 0  
info@ols-bio.de • www.ols-bio.com

OLS OMNI Life Science GmbH  
Laufenstraße 90  
4053 Basel, Switzerland  
Phone: +49 421 27 61 69 0  
info@ols-bio.ch • www.ols-bio.com

