

**CASY<sup>VIVO</sup> Cell Counter & Analyzer - Powering Breakthroughs in Cell Research**

**High-Precision Growth Analysis of Pterin-Auxotrophic Trypanosomatids**

Ong, H. B. et al. (2015). Comparative metabolism of conjugated and unconjugated pterins in Crithidia, Leishmania and African trypanosomes. *PLoS Neglected Tropical Diseases*, 9(12), e0004245. DOI: 10.1371/journal.pntd.0004245.

<i>Leishmania major</i> ; <i>Crithidia fasciculata</i> ; <i>Trypanosome brucei</i>	
Index	BP1
Standardization	
Counting	X
Viability	X
Volume	

**The Challenge:** Investigating why kinetoplastid parasites are pterin auxotrophs and determining the specific structural requirements of pterin analogues needed to support their growth and survival.

**CASY's Contribution:** CASY enabled precise daily monitoring of cell densities to determine half-maximal growth GC<sub>50</sub> values for various pterins. *T. brucei* and *C. fasciculata* have nearly identical growth preferences.


**Key Benefits to Researchers:**

**Quantitative Precision:** CASY provided the exact automated counts required to calculate precise half-maximal growth GC<sub>50</sub> values for diverse pterin analogues.

**Label-Free Efficiency:** The system allowed for continuous, dye-free tracking of growth kinetics over 14-day periods, avoiding interference with delicate metabolic experiments.

**Survival Validation:** Automated monitoring definitively confirmed the timing of growth cessation and viability loss in the absence of essential pterin nutrients.

**Figure 4 (Cumulative growth curves).** This graph displays daily cell densities measured by the CASY system, illustrating the halt in growth and progressive loss of viability when parasites are deprived of biopterin.



CASY was **essential** for quantifying the nutritional requirements of these parasites. Its ability to provide precise, automated daily counts allowed researchers to definitively map the metabolic dependency of *T. brucei* on biopterin, validating pteridine reductase 1 as a key survival factor and a potential drug target.